

Donkey Anti-Mouse IgG H&L (PE) preadsorbed ab7003

[9 References](#) [画像数 1](#)

製品の概要

製品名	Donkey Anti-Mouse IgG H&L (PE) preadsorbed
由来種	Donkey
ターゲット生物種	Mouse
アプリケーション	適用あり: IHC-P, IHC-Fr, Immunomicroscopy, Flow Cyt, ICC/IF
吸着処理血清	Chicken, Cow, Goat, Guinea pig, Hamster, Horse, Human, Rabbit, Rat, Sheep more details
免疫原	Mouse IgG whole molecule
標識	PE. Ex: 488nm, Em: 575nm

製品の特性


製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA
精製度	Immunogen affinity purified
特記事項 (精製)	Prepared from monospecific antiserum by immunoaffinity chromatography using Mouse IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation
特記事項 (標識)	R-Phycoerythrin (R-PE) Molecular Weight 240,000 daltons) Absorption Wavelength: 488 nm Emission Wavelength: 575 nm OD _{566nm} /OD _{280nm} : 2.8
ポリモノ	ポリクローナル
アイソタイプ	IgG
特記事項	Can be used in immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity.

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab7003の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		Use at an assay dependent dilution.
IHC-Fr		Use at an assay dependent dilution.
Immunomicroscopy		Use at an assay dependent dilution.
Flow Cyt		1/100 - 1/250.
ICC/IF		1/100 - 1/250.

画像



ab7003 was used with the primary antibody **ab33315** in Flow Cyt.
See Abreview on **ab33315**.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors