

Donkey Anti-Rabbit IgG H&L (10nm Gold) ab39597

2 References

製品の概要

製品名	Donkey Anti-Rabbit IgG H&L (10nm Gold)
由来種	Donkey
ターゲット生物種	Rabbit
アプリケーション	適用あり: Electron Microscopy
免疫原	Rabbit IgG whole molecule.
標識	Gold 10nm

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.097% Sodium azide Constituents: PBS, 1% BSA
精製度	Immunogen affinity purified
特記事項(精製)	This product has been immuno-affinity purified and immuno cross-absorbed to reduce non-specific reactions. The activity of each lot is determined using a dot-spot test system as described by Moeremans et al., J. Immunol. Methods, 74, (1984), 353.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG
特記事項	This product is an immunogold reagent built around colloidal gold particles of 10nm. The particle population is monodisperse and thus shows minimal size variation and overlap. Typically, the coefficient of variance for the 10nm particle size conjugate is less than 10 %.

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab39597の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Electron Microscopy		

追加情報

Electron Microscopy: Use at an assay dependent dilution.
Light Microscopy: Use at an assay dependent dilution.

Not yet tested in other applications.
Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors