

Donkey Anti-Rabbit IgG H&L (12nm Gold) preadsorbed ab105295

1 References

製品の概要

製品名	Donkey Anti-Rabbit IgG H&L (12nm Gold) preadsorbed
由来種	Donkey
ターゲット生物種	Rabbit
特異性	Recognizes rabbit IgG, H&L. X-Adsorbed and exhibits minimal crossreactivity to bovine, chicken, goat, guinea pig, Syrian hamster, equine, human, mouse, rat and sheep serum proteins.
アプリケーション	適用あり: Electron Microscopy, ELISA, WB
吸着処理血清	Chicken, Cow, Goat, Guinea pig, Horse, Human, Mouse, Rat, Sheep, Syrian hamster
標識	Gold 12nm

[more details](#)

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 8.50 Preservative: 0.05% Sodium azide Constituents: 0.164% Sodium phosphate, 0.381% Sodium borate, 1.5% BSA, 0.87% Sodium chloride
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab105295の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Electron Microscopy		1/20 - 1/40.
ELISA		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors