

Donkey Anti-Rat IgG H&L (DyLight® 650) preadsorbed ab102263

7 References

製品の概要

製品名	Donkey Anti-Rat IgG H&L (DyLight® 650) preadsorbed
由来種	Donkey
ターゲット生物種	Rat
特異性	By immunoelectrophoresis and ELISA this antibody reacts specifically with rat IgG and with light chains common to other rat immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. Less than 1% cross reactivity to bovine, chicken, goat, horse, human, mouse, rabbit and sheep IgG was detected. This antibody may cross react with IgG from other species.
アプリケーション	適用あり: Flow Cyt, ICC, IHC-P
吸着処理血清	Chicken, Cow, Goat, Horse, Human, Mouse, Rabbit, Sheep more details
標識	DyLight® 650. Ex: 654nm, Em: 673nm

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: PBS, 0.2% BSA
精製度	Immunogen affinity purified
特記事項(精製)	Antiserum was cross adsorbed using bovine, chicken, goat, horse, human, mouse, rabbit and sheep immunosorbents to remove cross reactive antibodies. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 650.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG
特記事項	DyLight® 650 replaces DyLight® 649.

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab102263の使用に適用されず**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		1/50 - 1/200.
ICC		1/50 - 1/500.
IHC-P		1/50 - 1/500. DyLight® 650 replaces DyLight® 649. The DyLight® 650 excitation is 652nm and the emission is 672nm. DyLight® 650 hence provides the same far-red fluorescence and photostability as the DyLight® 649 dye. Please refer to the vial in order to check whether the product you received is DyLight® 649 or DyLight® 650. This information however will have no impact on your

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors