

### Donkey Anti-Sheep IgG H&L (HRP) ab97125

★★★★☆ **1 Abreviews** **10 References**

#### 製品の概要

製品名	Donkey Anti-Sheep IgG H&L (HRP)
由来種	Donkey
ターゲット生物種	Sheep
特異性	By immunoelectrophoresis and ELISA this antibody reacts specifically with Sheep IgG and with light chains common to other Sheep immunoglobulins. No antibody was detected against non immunoglobulin serum proteins.
アプリケーション	<b>適用あり:</b> ICC, IHC-P, ELISA, WB
免疫原	Other Immunogen Type corresponding to Sheep IgG.
標識	HRP

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.8 Constituents: 0.2% BSA, 0.05% CMIT/MIT based preservative
精製度	Immunogen affinity purified
特記事項 (精製)	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Horse Radish Peroxidase (HRP).
特記事項 (標識)	Molar enzyme/ antibody protein ratio is 4:1
ポリモノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab97125の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC		Use at an assay dependent dilution.
IHC-P		1/200 - 1/5000.
ELISA		1/10000 - 1/100000. (primary)
WB	★★★★★ (1)	1/5000 - 1/50000. Colorimetric: 1/500 - 1/30000, Chemiluminescent: 1/10000 - 1/50000

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors