

# Goat Anti-Bird IgY H&L (HRP) ab112773

## 2 References

### 製品の概要

|          |   |
|----------|---|
| 製品名      | Goat Anti-Bird IgY H&L (HRP)  |
| 由来種      | Goat  |
| ターゲット生物種 | Bird  |
| 特異性      | By immunoelectrophoresis and ELISA ab112773 reacts with Dove, Duck, Sparrow, Chicken, Turkey, Emu, Ostrich, Quail, Macaw and Cockateil IgG. No antibody was detected against nonimmunoglobulin serum proteins. May cross react with IgG from other species. |
| アプリケーション | <b>適用あり:</b> WB, ELISA, IHC-P, ICC/IF   |
| 免疫原      | Other Immunogen Type corresponding to IgY.  |
| 標識       | HRP   |

### 製品の特性

|          |  |
|----------|--|
| 製品の状態    | Liquid   |
| 保存方法     | Shipped at 4°C. Store at +4°C.   |
| バッファー    | pH: 6.8<br>Constituents: 0.2% BSA, 0.05% CMIT/MIT based preservative   |
| 精製度      | Immunogen affinity purified  |
| 特記事項(精製) | Goats were immunized with purified IgG from Dove, Duck, Sparrow and Chicken serum. ab112773 was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to horseradish peroxidase (HRP). |
| 特記事項(標識) | Molar enzyme/antibody protein ratio is 4:1.  |
| ポリ/モノ    | ポリクローナル  |
| アイソタイプ   | IgG  |

### アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab112773の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項   |
|----------|-----------|--|
| WB       |           | 1/1000 - 1/10000. colorimetric detection; 1/1000 - 1/30000 chemiluminescent detection. |
| ELISA    |           | 1/10000 - 1/100000. as a primary.  |
| IHC-P    |           | 1/200 - 1/500.   |
| ICC/IF   |           | 1/200 - 1/500.   |

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors