

Goat Anti-Chicken IgY H&L (25nm Gold) ab41494

製品の概要

製品名	Goat Anti-Chicken IgY H&L (25nm Gold)
由来種	Goat
ターゲット生物種	Chicken
アプリケーション	<b>適用あり:</b> Electron Microscopy
免疫原	Chicken IgG whole molecule.
標識	Gold 25nm

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.097% Sodium azide Constituents: PBS, 1% BSA
精製度	Immunogen affinity purified
特記事項 (精製)	This product has been immuno-affinity purified and immuno cross-absorbed to reduce non-specific reactions. The activity of each lot is determined using a dot-spot test system [Moeremans et al., J. Immunol., Methods 74, (1984), 353.]
特記事項 (標識)	This product is an immunogold reagent built around colloidal gold particles of 6nm. The particle population is monodisperse and thus shows minimal size variation and overlap. Typically, the coefficient of variance for the 6nm particle size conjugate is less than 12%.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee      Abpromise保証は、次のテスト済みアプリケーションにおけるab41494の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Electron Microscopy		

## 追加情報

Electron Microscopy: Use at an assay dependent dilution. Electron Microscopy: Use at an assay dependent dilution. The method of production of this product allows the small particle size of 6nm to be used instead of 5nm with the same sensitivity but increased visibility in the standard transmission electron microscope (See also van Bergen en Henegouwen and Leunissen, Histochemistry, 85,(1986), 81.)

Light Microscopy: Use at an assay dependent dilution.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors