

Goat Anti-Chicken IgY H&L (FITC) ab46969

6 References 画像数 1

製品の概要

製品名	Goat Anti-Chicken IgY H&L (FITC)
由来種	Goat
ターゲット生物種	Chicken
アプリケーション	適用あり: ICC/IF, IHC-Fr
免疫原	Chicken IgY
標識	FITC. Ex: 493nm, Em: 528nm

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7 Preservative: 0.02% Sodium azide Constituents: PBS, 0.2% BSA Isotonic 0.9%, w/v)
精製度	Affinity purified
特記事項 (精製)	Ab46969 is specific to chicken IgY and will not cross react with mammalian IgG.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、**次のテスト済みアプリケーションにおけるab46969の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		

アプリケーション	Abreviews	特記事項
IHC-Fr		

追加情報


ICC/IF: 1/250.

IHC-Fr: Use at an assay dependent dilution (PMID 18583367).

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

画像



ab4674 staining GFAP in human brain cells by Flow Cytometry.

Cells were cultured in stem cell media RHB-A and collected using

0.05% Trypsin. The sample was incubated with the primary

antibody (1/100 in PBS + 10% FCS) for 20 minutes at 25°C. A

FITC-conjugated Goat anti-chicken IgY (ab46969) (1/100) was used as the secondary antibody.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors