

## Product datasheet

# Goat F(ab')<sub>2</sub> Anti-Human IgG - Fc (DyLight® 488), pre-adsorbed ab98590

画像数 1

### 製品の概要

製品名	Goat F(ab') <sub>2</sub> Anti-Human IgG - Fc (DyLight® 488), pre-adsorbed
由来種	Goat
ターゲット生物種	Human
アプリケーション	適用あり: Flow Cyt, ICC/IF, IHC-P
吸着処理血清	Mouse, Rat
標識	DyLight® 488. Ex: 493nm, Em: 518nm

[more details](#)

### 製品の特性


製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: PBS, 0.2% BSA
精製度	Immunogen affinity purified
特記事項(精製)	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. F(ab') <sub>2</sub> fragment were generated using a pepsin digestion. Fc fragments and whole IgG molecules have been removed. Fragments were conjugated to DyLight® 488.
ポリモノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

**The Abpromise guarantee** **Abpromise保証は、次のテスト済みアプリケーションにおけるab98590の使用に適用されず**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		1/50 - 1/200.
ICC/IF		1/50 - 1/500.
IHC-P		1/50 - 1/500.

## 画像



Goat F(ab')<sub>2</sub> Anti-Human IgG - Fc (DyLight® 488),  
pre-adsorbed (ab98590)

Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors