

Goat F(ab')<sub>2</sub> Anti-Mouse IgG - (Fab)'<sub>2</sub> (PE), pre-adsorbed ab5889

2 References

製品の概要	
製品名	Goat F(ab') <sub>2</sub> Anti-Mouse IgG - (Fab)' <sub>2</sub> (PE), pre-adsorbed
由来種	Goat
ターゲット生物種	Mouse
アプリケーション	適用あり: ICC/IF, Immunomicroscopy, Flow Cyt
吸着処理血清	Cow, Horse, Human, Rabbit, Rat, Sheep
免疫原	Mouse IgG F(ab) <sub>2</sub> fragment
標識	PE. Ex: 488nm, Em: 575nm

[more details](#)

製品の特性	
製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.878% Sodium chloride, 1% BSA, 0.424% Tripotassium orthophosphate
特記事項 (精製)	BSA is IgG and protease free This product was prepared from monospecific antiserum by immunoaffinity chromatography using Mouse IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation.
特記事項 (標識)	R-Phycoerythrin (R-PE) (Molecular Weight 240,000 daltons) from seaweed Absorption Wavelength: 490 nm, 545 nm and 565 nm Emission Wavelength: 580 nm OD <sub>566nm</sub> /OD <sub>280nm</sub> : 2.6 Phycobiliprotein Concentration: 0.5 mg/ml (determined by absorbance = 82.0 at 565 nm for a 1% solution for only those R-PE molecules to which at least one molecule of active antibody is bound).
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

**The Abpromise guarantee** Abpromise保証は、 次のテスト済みアプリケーションにおけるab5889の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		1/100 - 1/250.
Immunomicroscopy		Use at an assay dependent dilution. Suitable for other antibody based fluorescent assays requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity.
Flow Cyt		1/100 - 1/250. The maximum amount of reagent required to stain $1 \times 10^6$ cells in flow cytometry is approximately 1.0 µg of antibody conjugate. Lesser amounts of reagent may be sufficient for staining.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors