

### Goat F(ab')<sub>2</sub> Anti-Mouse IgG - H&L (HRP), pre-adsorbed ab5870

#### 10 References

#### 製品の概要

製品名	Goat F(ab') <sub>2</sub> Anti-Mouse IgG - H&L (HRP), pre-adsorbed
由来種	Goat
ターゲット生物種	Mouse
アプリケーション	<b>適用あり:</b> ICC/IF, Dot blot, WB, ELISA, Electron Microscopy, IHC
吸着処理血清	Cow, Horse, Human, Rabbit, Rat, Sheep <a href="#">more details</a>
免疫原	Mouse IgG whole molecule.
標識	HRP

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.50 Preservative: 0.01% Gentamicin sulphate Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA  BSA is IgG and protease free
特記事項 (精製)	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Mouse IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation.
特記事項 (標識)	Horseradish Peroxidase (HRP) (Molecular Weight 40,000 daltons)
ポリモノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab5870の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		1/500 - 1/2500.
Dot blot		Use at an assay dependent dilution.
WB		1/1000 - 1/10000.
ELISA		1/60000.
Electron Microscopy		Use at an assay dependent dilution.
IHC		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors