

Goat F(ab')₂ Anti-Mouse IgG+IgM+IgA H&L (APC) preadsorbed
ab130775

| | |
|----------|--|
| 製品の概要 | |
| 製品名 | Goat F(ab') ₂ Anti-Mouse IgG+IgM+IgA H&L (APC) preadsorbed |
| 由来種 | Goat |
| ターゲット生物種 | Mouse |
| 特異性 | ab130775 reacts with the heavy and light chains of mouse IgM, IgG1, IgG2a, IgG2b IgG3, and IgA. Goat F(ab') ₂ Anti-Mouse IgG+IgM+IgA H&L (Allophycocyanin) preadsorbed shows minimal cross reactivity with Human immunoglobulins. |
| アプリケーション | 適用あり: ICC/IF, ELISA, WB, Dot blot, IHC-P, IHC-Fr |
| 吸着処理血清 | Human more details |
| 免疫原 | Full length native protein (purified) corresponding to Mouse IgG + IgM + IgA. Mouse Ig |
| 標識 | APC. Ex: 645nm, Em: 660nm |

| | |
|-----------|--|
| 製品の特性 | |
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C. Store In the Dark. |
| バッファー | pH: 7.3 Preservative: 0.1% Sodium azide Constituent: PBS Sucrose is added as a stabilizing agent. |
| 特記事項 (精製) | ab130775 was purified by gel filtration chromatography of pepsin digested antibody. |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |
| 特記事項 | ab130775 is adsorbed against pooled Human sera and purified Human paraproteins. |

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab130775の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|--|
| ICC/IF | | Use 0.1µg for 10 ⁶ cells. |
| ELISA | | Use at an assay dependent dilution. |
| WB | | Use at an assay dependent dilution. |
| Dot blot | | Use at an assay dependent dilution. |
| IHC-P | | Use at an assay dependent dilution. |
| IHC-Fr | | Use at an assay dependent concentration. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors