

Goat F(ab')₂ Anti-Mouse IgG+IgM+IgA - H&L (Biotin), pre-adsorbed ab6005

1 References

製品の概要

製品名	Goat F(ab') ₂ Anti-Mouse IgG+IgM+IgA - H&L (Biotin), pre-adsorbed
由来種	Goat
ターゲット生物種	Mouse
特異性	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-biotin and anti-Goat Serum. No reaction was observed against Human Serum proteins. This product is suitable for the detection of all Mouse immunoglobulin classes, isotypes and chain combinations.
アプリケーション	適用あり: ICC/IF, Dot blot, WB, ELISA, IHC, Immunomicroscopy
吸着処理血清	Human more details
免疫原	Mouse IgG, IgA and IgM whole molecules
標識	Biotin

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA BSA is IgG and protease free
特記事項 (精製)	This product was prepared from polyspecific antiserum by immunoaffinity chromatography using antigens coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
特記事項 (標識)	Biotinamidocaproate N-Hydroxysuccinimide Ester (BAC) Biotin/Protein Ratio: 10-20 BAC molecules per Goat IgG molecule
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab6005の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		1/2000 - 1/20000.
Dot blot		Use at an assay dependent dilution.
WB		1/5000 - 1/20000.
ELISA		1/1.2e+006.
IHC		Use at an assay dependent dilution.
Immunomicroscopy		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors