

Goat F(ab')₂ Anti-Mouse IgM mu chain (TRITC) ab5928

2 References

製品の概要

| | |
|----------|---|
| 製品名 | Goat F(ab') ₂ Anti-Mouse IgM mu chain (TRITC) |
| 由来種 | Goat |
| ターゲット生物種 | Mouse |
| 特異性 | Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum, Mouse IgM and Mouse Serum. No reaction was observed against anti-Pepsin or anti-Goat IgG F(c). Specificity was confirmed by ELISA at less than 1% cross reactivity against other mouse heavy or light chain isotypes. |
| アプリケーション | 適用あり: Functional Studies, ICC/IF, Immunomicroscopy, Flow Cyt |
| 免疫原 | Mouse IgM mu chain |
| 標識 | TRITC. Ex: 547nm, Em: 572nm |

製品の特性

| | |
|-----------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| バッファー | pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA BSA is IgG and protease free |
| 特記事項 (精製) | This product was prepared from monospecific antiserum by immunoaffinity chromatography using Mouse IgM coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation. |
| 特記事項 (標識) | Label: Tetramethylrhodamine isothiocyanate (TRITC) (Molecular Weight 444 daltons) Absorption Wavelength: 550 nm Emission Wavelength: 570 nm Fluorochrome/Protein Ratio: 2.3 moles TRITC per mole of Goat IgG F(ab') ₂ |
| ポリモノ | ポリクローナル |
| アイソタイプ | IgG |

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab5928の使用に適用されます**
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

| アプリケーション | Abreviews | 特記事項 |
|--------------------|-----------|-------------------------------------|
| Functional Studies | | 1/10000 - 1/50000. |
| ICC/IF | | 1/1000 - 1/5000. |
| Immunomicroscopy | | Use at an assay dependent dilution. |
| Flow Cyt | | 1/500 - 1/2500. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors