

Goat F(ab')₂ Anti-Rat IgG - (Fab')₂ (AP), pre-adsorbed ab6518

製品の概要

製品名	Goat F(ab') ₂ Anti-Rat IgG - (Fab') ₂ (AP), pre-adsorbed
由来種	Goat
ターゲット生物種	Rat
特異性	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Alkaline Phosphatase, anti-Goat Serum, Rat IgG F(ab) ₂ and Rat Serum. No reaction was observed against anti-Pepsin, anti-Goat IgG F(c), Rat Fc or Bovine, Horse, and Human Serum Proteins
アプリケーション	適用あり: ELISA, IHC, ICC/IF, Dot blot, WB
吸着処理血清	Cow, Horse, Human more details
免疫原	Rat IgG F(ab) ₂ fragment
標識	Alkaline Phosphatase

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.05% Sodium azide Constituents: 0.79% Tris HCl, 0.88% Sodium chloride, 1.5% BSA, 0.01% Magnesium chloride, 0.0014% Zinc chloride, 50% Glycerol (glycerin, glycerine) BSA is IgG and protease free
特記事項 (精製)	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rat IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation.
特記事項 (標識)	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、**次のテスト済みアプリケーションにおけるab6518の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
ELISA		1/5000 - 1/50000.
IHC		Use at an assay dependent dilution.
ICC/IF		1/200 - 1/1000.
Dot blot		Use at an assay dependent dilution.
WB		1/5000 - 1/50000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors