

Goat Anti-Human IgA (SureLight® APC) ab99897

1 References [画像数 1](#)

製品の概要

製品名	Goat Anti-Human IgA (SureLight® APC)
由来種	Goat
ターゲット生物種	Human
アプリケーション	適用あり: Flow Cyt
免疫原	Full length native protein (purified) corresponding to IgA.
標識	SureLight® APC. Ex: 652nm, Em: 657nm

製品の特性


製品の状態	Lyophilized:Reconstitute with 1 ml of ddH2O
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 8.20 Constituents: 0.1% BSA, 0.87% Sodium chloride, 0.121% Tris, 1.71% Sucrose
精製度	Immunogen affinity purified
特記事項(標識)	Fluor:Protein ~2:1 Visible absorption maxima 652 nm Emission maximum 657 nm
ポリ/モノ	ポリクローナル
アイソタイプ	IgG
特記事項	After product has been reconstituted, product should be stored at 2-8°C in the dark and be used within 3 months. If further dilution of the conjugate is required, use diluted material within one week.

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab99897の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		Use a concentration of 3 µg/ml.

画像



Flow Cytometry - Goat Anti-Human IgA (SureLight®
APC) (ab99897)

Fluorescence excitation and emission spectra of Allophycocyanin in 100mM sodium phosphate (pH7.2) + 1mM EDTA and 1mM sodium azide. The emission scan was taken with excitation at 630nm. The excitation was taken with an emission at 660nm. The curves were normalized to equalize peak heights.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors