

Goat Anti-Human IgG (Fab')₂ (HRP) ab87422

6 References [画像数 1](#)

医薬用外劇物

製品の概要

製品名	Goat Anti-Human IgG (Fab') ₂ (HRP)
由来種	Goat
ターゲット生物種	Human
特異性	ab87422 is specific for the Fab and F(ab') ₂ portion of human IgG by IEP. It reacts with other human immunoglobulins through common light chain reactivity.
アプリケーション	適用あり: ELISA, WB
免疫原	Other Immunogen Type corresponding to IgG. Human serum.
標識	HRP

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 0.1% BSA, PBS
精製度	Protein G purified
特記事項(精製)	This antibody has been cross adsorbed against the Fc domain and will react with less than 1 % of the Fc domain of human IgG.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG


アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab87422の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		1/1000 - 1/10000.

アプリケーション	Abreviews	特記事項
WB		1/1000 - 1/10000.

画像



Characterisation of virus activity and functional expression of monoclonal antibodies.

(Panel D) IgG1 antibody expression, assessed by anti-IgG1 ELISA, in cellular supernatants of carcinoma cell lines infected for 24-96hrs with NG-135 at 1ppc.

Media was removed from infected or control cells, clarified by centrifugation and diluted 1 in 2 in 3% BSA/PBS for storage at -20°C. For quantification of human IgG1, 96 well plates were coated overnight, 4°C, with mouse monoclonal anti-human IgG1 Fc antibody (**ab1927**, Abcam) at a 1:1000 dilution in carbonate/bicarbonate buffer. Plates were blocked in 3% BSA/PBS, washed and samples and standards added to the plate for 1 hr and RT. Secondary detection was carried out using a HRP conjugated goat anti-human IgG Fab (ab87422, Abcam) incubated for 1 hr, RT.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors