

Goat Anti-Human IgG H&L (HRP) ab81202

4 References

製品の概要

製品名	Goat Anti-Human IgG H&L (HRP)
由来種	Goat
ターゲット生物種	Human
特異性	No crossreaction with human IgM and IgA in ELISA No crossreaction with albumin fraction of human serum. No crossreaction with albumin fraction of other animal species.
アプリケーション	適用あり: WB, ELISA
免疫原	Human IgG.
標識	HRP

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.50 Constituents: 50% Glycerol, 0.58% Sodium chloride, 0.164% Sodium phosphate
精製度	Immunogen affinity purified
特記事項 (精製)	From the crude polyclonal the crossreactive antibodies were extracted by incubation with Sepharose bound human IgA and IgM. Specific antibodies were absorbed by incubation with Sepharose bound human IgG. Specific antibodies were eluted by acidic buffer at pH 2.5 followed by neutralization and dialysis. After repeated binding with immobilized human IgG a minimum of 65% protein bound.
特記事項 (標識)	The purified polyclonal was conjugated to horseradish peroxidase according to the periodate method followed by gel filtration and ion exchange chromatography to clear unbound conjugate.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab81202の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/100.
ELISA		1/1000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors