

### Goat Anti-Human IgG H&L preadsorbed ab7148

#### 2 References

#### 製品の概要

製品名	Goat Anti-Human IgG H&L preadsorbed
由来種	Goat
ターゲット生物種	Human
アプリケーション	<b>適用あり:</b> WB, IHC-P, IHC-Fr, IP, ICC/IF, Immunodiffusion, Conjugation, ELISA
吸着処理血清	Chicken, Cow, Goat, Guinea pig, Hamster, Horse, Mouse, Rabbit, Rat, Sheep <a href="#">more details</a>
免疫原	Human IgG whole molecule

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.01% Sodium azide Constituents: 0.88% Sodium chloride, 0.42% Potassium phosphate
精製度	Affinity purified
特記事項 (精製)	This product was prepared from normal serum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis.
ポリ/モノ	ポリクローナル

#### アプリケーション

**The Abpromise guarantee** **Abpromise保証は、次のテスト済みアプリケーションにおけるab7148の使用に適用されます**  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/2000 - 1/10000.
IHC-P		Use at an assay dependent concentration.

アプリケーション	Abreviews	特記事項
IHC-Fr		Use at an assay dependent concentration.
IP		Use at an assay dependent concentration.
ICC/IF		1/1000 - 1/5000.
Immunodiffusion		Use at an assay dependent concentration.
Conjugation		Use at an assay dependent concentration.
ELISA		1/25000 - 1/120000.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors