

Goat Anti-Human IgG+IgM+IgA H&L ab102416

1 References

製品の概要

製品名	Goat Anti-Human IgG+IgM+IgA H&L
由来種	Goat
ターゲット生物種	Human
特異性	By immunoelectrophoresis and ELISA this antibody reacts with human IgG, IgA and IgM. No antibody was detected against non-immunoglobulin serum proteins. This antibody may cross react with immunoglobulins from other species.
アプリケーション	適用あり: ICC, IHC-P, ELISA, WB
免疫原	Other Immunogen Type corresponding to IgG + IgM + IgA.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 6.8 Preservative: 0.09% Sodium azide Constituent: PBS
精製度	Immunogen affinity purified
特記事項(精製)	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads.
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab102416の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC		Use at an assay dependent dilution.

アプリケーション	Abreviews	特記事項
IHC-P		1/200 - 1/2000.
ELISA		1/100 - 1/30000. ELISA - Coating 1/100 - 1/500; ELISA - Primary 1/1000 - 1/30000
WB		1/1000 - 1/30000. WB - Colorimetric 1/1000 - 1/10000; WB - Chemiluminescent 1/1000 - 1/30000

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors