

Goat Anti-Human IgG+IgM+IgA H&L (Biotin) ab102418

1 References

製品の概要

製品名	Goat Anti-Human IgG+IgM+IgA H&L (Biotin)
由来種	Goat
ターゲット生物種	Human
特異性	By immunoelectrophoresis and ELISA this antibody reacts with human IgG, IgA and IgM. No antibody was detected against non-immunoglobulin serum proteins. This antibody may cross react with immunoglobulins from other species.
アプリケーション	適用あり: IHC-P, ICC/IF, ELISA, WB
免疫原	Other Immunogen Type corresponding to IgG + IgM + IgA.
標識	Biotin

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: PBS, 0.2% BSA
精製度	Immunogen affinity purified
特記事項(精製)	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Biotin.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、**次のテスト済みアプリケーションにおけるab102418の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/500 - 1/5000.
ICC/IF		1/200 - 1/1000.
ELISA		1/10000 - 1/200000. ELISA - Primary
WB		1/10000 - 1/200000. WB - Chemiluminescent

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors