

Goat Anti-Human IgM H&L (FITC) ab8497

製品の概要

製品名	Goat Anti-Human IgM H&L (FITC)
由来種	Goat
ターゲット生物種	Human
アプリケーション	適用あり: ICC/IF
標識	FITC. Ex: 493nm, Em: 528nm

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.50 Constituents: 0.164% Sodium phosphate, 50% Glycerol, 0.58% Sodium chloride
精製度	Immunogen affinity purified
特記事項 (精製)	From the crude polyclonal the crossreactive antibodies were extracted by incubation with Sepharose-bound human IgA and IgG. Specific antibodies were absorbed by incubation with Sepharose-bound human IgM. Specific antibodies were eluted by acidic buffer at pH 2.5 followed by neutralisation and dialysis. After repeated binding with immobilized human IgM a minimum of 65% protein bound.
特記事項 (標識)	The purified polyclonal was conjugated to fluorescein isothiocyanate isomer I followed by gel-filtration and ion-exchange chromatography to clear unbound conjugate.
ポリモノ	ポリクローナル
軽鎖の種類	unknown
特記事項	F/P ratio is 4-6 mol Fluorescein per 1 mol of goat IgM. There is low background fluorescence in normal tissue in contrast to cases of auto immune disease where IgM is deposited.

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab8497の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		

追加情報

For immunohistochemical staining of human tissue sections a 1:40 dilution is recommended.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors