

Goat Anti-Monkey IgG H&L ab112764

製品の概要

製品名	Goat Anti-Monkey IgG H&L
由来種	Goat
ターゲット生物種	Monkey
特異性	By immunoelectrophoresis and ELISA ab112764 reacts specifically with monkey IgG and with light chains common to other monkey immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. ab112764 may cross react with IgG from other species.
アプリケーション	適用あり: WB, ELISA, IHC-P, ICC/IF
免疫原	Other Immunogen Type corresponding to IgG.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 6.8 Preservative: 0.1% Sodium azide Constituent: 99% PBS
精製度	Immunogen affinity purified
特記事項 (精製)	ab112764 was isolated by affinity chromatography using antigen coupled to agarose beads.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab112764の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/1000 - 1/10000. colorimetric detection; 1/1000 - 1/30000 chemiluminescent detection.

アプリケーション	Abreviews	特記事項
ELISA		1/100 - 1/500. for coating; 1/1000 - 1/30000 as a primary.
IHC-P		1/200 - 1/2000.
ICC/IF		1/200 - 1/2000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors