

# Goat Anti-Monkey IgG H&L (FITC) ab112766

## 2 References

### 製品の概要

製品名	Goat Anti-Monkey IgG H&L (FITC)
由来種	Goat
ターゲット生物種	Monkey
特異性	By immunoelectrophoresis and ELISA ab112766 reacts specifically with monkey IgG and with light chains common to other monkey immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. ab112766 may cross react with IgG from other species.
免疫原	Other Immunogen Type corresponding to IgG.
標識	FITC. Ex: 493nm, Em: 528nm

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, 99% PBS
精製度	Immunogen affinity purified
特記事項 (精製)	ab112766 was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to fluorescein isothiocyanate (FITC).
特記事項 (標識)	F/P ratio is 4.2.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors