

### Goat Anti-Mouse IgG Fc (DyLight® 650) ab97018

[8 References](#) [画像数 1](#)

#### 製品の概要

|          |   |
|----------|---|
| 製品名      | Goat Anti-Mouse IgG Fc (DyLight® 650)   |
| 由来種      | Goat  |
| ターゲット生物種 | Mouse   |
| 特異性      | By immunoelectrophoresis and ELISA this antibody reacts specifically with Mouse IgG-Fc. Cross reactivity with other immunoglobulins and light chains is less than 0.1%. |
| アプリケーション | <b>適用あり:</b> IHC-P, ICC/IF, Flow Cyt  |
| 標識       | DyLight® 650. Ex: 654nm, Em: 673nm  |

#### 製品の特性


|           |  |
|-----------|--|
| 製品の状態     | Liquid   |
| 保存方法      | Shipped at 4°C. Store at +4°C.   |
| バッファー     | pH: 6.8<br>Preservative: 0.09% Sodium azide<br>Constituents: 0.2% BSA, PBS   |
| 精製度       | Immunogen affinity purified  |
| 特記事項 (精製) | Antiserum was solid phase adsorbed to ensure class specificity. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 650. |
| ポリ/モノ     | ポリクローナル  |
| アイソタイプ    | IgG  |
| 特記事項      | DyLight® 650 replaces DyLight® 649.  |

#### アプリケーション

**The Abpromise guarantee** Abpromise保証は、次のテスト済みアプリケーションにおけるab97018の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項          |
|----------|-----------|---------------|
| IHC-P    |           | 1/50 - 1/500. |
| ICC/IF   |           | 1/50 - 1/500. |
| Flow Cyt |           | 1/50 - 1/200. |

## 画像



Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors