

### Goat Anti-Mouse IgG Fc (FITC) ab97264

5 References [画像数 1](#)

#### 製品の概要

製品名	Goat Anti-Mouse IgG Fc (FITC)
由来種	Goat
ターゲット生物種	Mouse
特異性	By immunoelectrophoresis and ELISA this antibody reacts specifically with Mouse IgG-Fc. Cross reactivity with other immunoglobulins and light chains is less than 0.1%.
アプリケーション	<b>適用あり:</b> IHC-P, ICC/IF
標識	FITC. Ex: 493nm, Em: 528nm

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, PBS
精製度	Immunogen affinity purified
特記事項 (精製)	Antiserum was solid phase adsorbed to ensure class specificity. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to FITC.
特記事項 (標識)	F/P ratio is 4.4
ポリモノ	ポリクローナル
アイソタイプ	IgG


#### アプリケーション

**The Abpromise guarantee** **Abpromise保証は、次のテスト済みアプリケーションにおけるab97264の使用に適用されます**  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/50 - 1/500.

アプリケーション	Abreviews	特記事項
ICC/IF		1/50 - 1/500.

## 画像



A431 cell nucleus (DAPI signal) showing three gamma-H2AX foci structures.

Primary antibody H2A.X (phospho S139) antibody [3F2], [ab22551](#), 100ug.

Secondary antibody Mouse IgG-Fc (FITC), ab97264, 1mg.

Immunocytochemistry/ Immunofluorescence - Goat

Anti-Mouse IgG Fc (FITC) (ab97264)

This image is courtesy of Jorge E. Gonzalez (CPHR, La Habana, Cuba), Joan F. Barquinero (UAB, Barcelona, Spain) and Jessica Martinez (UAB, Universitat Autònoma de Barcelona, Spain).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors