

### Goat Anti-Mouse IgG H&L (APC) preadsorbed ab130782

#### 4 References

#### 製品の概要

製品名	Goat Anti-Mouse IgG H&L (APC) preadsorbed
由来種	Goat
ターゲット生物種	Mouse
特異性	ab130782 reacts with the heavy and light chains of mouse IgG1, IgG2a, IgG2b and IgG3 and with the light chains of mouse IgM and IgA as demonstrated by ELISA, FLISA, and flow cytometry. Goat Anti-Mouse IgG H&L (Allophycocyanin) preadsorbed exhibits minimal cross reactivity with Human immunoglobulins.
アプリケーション	<b>適用あり:</b> ICC/IF, ELISA, WB
吸着処理血清	Human <a href="#">more details</a>
免疫原	Mouse IgG paraproteins
標識	APC. Ex: 645nm, Em: 660nm

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Store In the Dark.
バッファー	pH: 7.3 Preservative: 0.1% Sodium azide Constituents: PBS, Sucrose
特記事項 (精製)	ab130782 was purified by affinity chromatography on pooled Mouse IgG covalently linked to agarose.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG
特記事項	ab130782 is cross absorbed against pooled Human sera and purified Human paraproteins.

#### アプリケーション

**The Abpromise guarantee** **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab130782の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		Use 0.1µg for 10 <sup>6</sup> cells.
ELISA		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors