

### Goat Anti-Mouse IgG H&L (HRP) preadsorbed ab97040

★★★★★ [2 Abreviews](#) [252 References](#)

#### 製品の概要

製品名	Goat Anti-Mouse IgG H&L (HRP) preadsorbed
由来種	Goat
ターゲット生物種	Mouse
特異性	By immunoelectrophoresis and ELISA this antibody reacts specifically with Mouse IgG and with light chains common to other Mouse immunoglobulins. No antibody was detected against non immunoglobulin serum proteins. Reduced cross-reactivity to bovine, chicken, goat, horse, human, pig, rabbit and rat IgG was detected.
アプリケーション	<b>適用あり:</b> IHC-P, ELISA, WB
吸着処理血清	Chicken, Cow, Goat, Horse, Human, Pig, Rabbit, Rat, Sheep <a href="#">more details</a>
免疫原	Other Immunogen Type corresponding to Mouse IgG.
標識	HRP

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.8 Constituents: 0.2% BSA, PBS, 0.00063% 2-methyl-4-isothiazolin-3-one, 0.00063% 5-chloro-2-methyl-4-isothiazolin-3-one
精製度	Immunogen affinity purified
特記事項(精製)	Antiserum was cross adsorbed using bovine, chicken, horse, human, pig, rabbit and rat immunosorbents to remove cross reactive antibodies. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Horse Radish Peroxidase (HRP).
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

**The Abpromise guarantee** Abpromise保証は、 次のテスト済みアプリケーションにおけるab97040の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/200 - 1/5000.
ELISA		1/10000 - 1/100000.
WB	★★★★★ (2)	1/5000 - 1/50000. For colormetric detection 1/5000 - 1/30000. For chemiluminescent detection 1/10000 - 1/50000.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors