



Product datasheet

Goat Anti-Mouse IgG2a heavy chain ab9163

1 References

製品の概要

製品名	Goat Anti-Mouse IgG2a heavy chain
由来種	Goat
ターゲット生物種	Mouse
特異性	The affinity purified anti IgG2a only reacts with mouse IgG2a by immunoelectrophoresis and double diffusion when tested at 1 mg/ml against all mouse immunoglobulin classes and subclasses at 1 mg/ml.
アプリケーション	<b>適用あり:</b> ELISA, Immunodiffusion, Immunoelectrophoresis, WB
免疫原	Goats were immunized with highly purified mouse IgG2a in Freund’s adjuvant

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.20
精製度	Immunogen affinity purified
特記事項 (精製)	The antiserum was solid phase absorbed to make it IgG2a heavy chain specific and the anti IgG2A subsequently immunoaffinity purified off an IgG2a immunoabsorbant.
ポリモノ	ポリクローナル

アプリケーション

The Abpromise guarantee      **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab9163の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.
Immunodiffusion		Use at an assay dependent concentration.

アプリケーション	Abreviews	特記事項
Immunoelectrophoresis		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors