

### Goat Anti-Mouse IgG+IgM H&L (Biotin) preadsorbed ab47844

[3 References](#) [画像数 1](#)

#### 製品の概要

製品名	Goat Anti-Mouse IgG+IgM H&L (Biotin) preadsorbed
由来種	Goat
ターゲット生物種	Mouse
特異性	This antibody reacts with Mouse IgG and IgM heavy and light chains. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum, anti-Biotin, Mouse IgG and Mouse IgM.
アプリケーション	<b>適用あり:</b> ELISA, WB, ICC/IF, IHC-P
吸着処理血清	Cow, Horse, Human <a href="#">more details</a>
免疫原	Full length IgG and IgM molecules (Mouse)
標識	Biotin

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA
精製度	Immunogen affinity purified
特記事項 (精製)	Prepared from polyspecific antiserum by immunoaffinity chromatography using antigens coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
特記事項 (標識)	10-20 BAC (Biotinamidocaproate N-Hydroxysuccinimide Ester) molecules per IgG molecule
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション


The Abpromise guarantee

**Abpromise保証は、次のテスト済みアプリケーションにおけるab47844の使用に適用されません**

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
ELISA		1/20000 - 1/100000.
WB		1/2000 - 1/10000.
ICC/IF		1/1000 - 1/5000.
IHC-P		Use at an assay dependent dilution.

## 画像



Ab47844 was used at dilution 1/1000 with the primary antibody [ab24575](#) in ICC/IF. See the review on [ab24575](#).

Immunocytochemistry/ Immunofluorescence - Goat  
Anti-Mouse IgG+IgM H&L (Biotin) preadsorbed  
(ab47844)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors