

Goat Anti-Mouse IgM mu chain (HRP) ab97230

48 References 画像数 1

製品の概要

製品名	Goat Anti-Mouse IgM mu chain (HRP)
由来種	Goat
ターゲット生物種	Mouse
特異性	By immunoelectrophoresis and ELISA this antibody reacts specifically with Mouse IgM. Cross reactivity with other immunoglobulins and light chains is less than 0.1%.
アプリケーション	適用あり: ELISA, WB, ICC, IHC-P
標識	HRP

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.8 Constituents: 0.2% BSA, 0.05% CMIT/MIT based preservative
精製度	Immunogen affinity purified
特記事項 (精製)	Antiserum was solid phase adsorbed to ensure class specificity. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Horse Radish Peroxidase (HRP).
特記事項 (標識)	Molar enzyme/ antibody protein ratio is 4:1
ポリ/モノ	ポリクローナル
アイソタイプ	IgG


アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab97230の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		1/10000 - 1/100000. (Primary ELISA)

アプリケーション	Abreviews	特記事項
WB		1/2000 - 1/25000. (Colorimetric: 1/2000 - 1/20000; Chemiluminescence: 1/5000 - 1/25000)
ICC		Use at an assay dependent dilution.
IHC-P		1/200 - 1/2000.

画像



IHC image of Hsp70 staining in human lung formalin fixed paraffin embedded tissue section*. The section was pre-treated using pressure cooker heat mediated antigen retrieval with sodium citrate buffer (pH6) for 30mins. The section was incubated with **ab5442**, 1/2000 dilution overnight at +4°C. An HRP-conjugated secondary (ab97230, 1/2000 dilution) was used for 1hr at room temperature. The section was counterstained with haematoxylin and mounted with DPX.

The inset negative control image is secondary-only at 1/500 dilution.

*Tissue obtained from the Human Research Tissue Bank, supported by the NIHR Cambridge Biomedical Research Centre

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors