

Goat Anti-Mouse Anti-Rabbit IgG/IgM H&L (HRP polymer) ab2891

21 References

製品の概要

製品名	Goat Anti-Mouse Anti-Rabbit IgG/IgM H&L (HRP polymer)
由来種	Goat
ターゲット生物種	Mouse
特異性	This antibody reacts with heavy and light chains on mouse and rabbit IgG. It reacts with all mouse and rabbit IgG subclasses and with mouse and rabbit IgM.
アプリケーション	適用あり: IHC-P, IHC-Fr
標識	HRP polymer

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Constituents: 0.04% Proclin 300, 0.5% Proclin 950, 99% Water
特記事項 (標識)	HRP conjugated to a polymer backbone, which is collectively conjugated to the antibody. The effect is multiple HRP molecules linked to one single antibody.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG
特記事項	This new and innovative HRP-polymerization technology provides a significant increase in staining sensitivity when compared to other conventional HRP-conjugated secondary antibodies. Avidin-biotin blocking procedures are not necessary when using this product. The overall staining procedure uses one less reagent (streptavidin peroxidase) and rinsing step, eliminates avidin-biotin blocking procedures and may reduce overall staining time.

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab2891の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		Use at an assay dependent dilution.
IHC-Fr		Use at an assay dependent dilution. This product is provided in a pre-diluted, ready-to-use format. Apply for 20 - 30 minutes.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors