

### Goat Anti-Pig IgG H&L (HRP) ab102135

#### 4 References

#### 製品の概要

|          |   |
|----------|---|
| 製品名      | Goat Anti-Pig IgG H&L (HRP)   |
| 由来種      | Goat  |
| ターゲット生物種 | Pig   |
| 特異性      | By immunoelectrophoresis and ELISA this antibody reacts specifically with pig IgG and with light chains common to other pig immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. This antibody may cross react with IgG from other species. |
| アプリケーション | <b>適用あり:</b> IHC-P, ICC, ELISA, WB  |
| 標識       | HRP   |

#### 製品の特性

|           |   |
|-----------|---|
| 製品の状態     | Liquid  |
| 保存方法      | Shipped at 4°C. Store at +4°C.  |
| バッファー     | pH: 6.8<br>Constituents: 0.2% BSA, 0.05% CMIT/MIT based preservative  |
| 精製度       | Immunogen affinity purified   |
| 特記事項 (精製) | This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to HRP. |
| ポリ/モノ     | ポリクローナル   |
| アイソタイプ    | IgG   |

#### アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab102135の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項            |
|----------|-----------|-----------------|
| IHC-P    |           | 1/200 - 1/5000. |

| アプリケーション | Abreviews | 特記事項  |
|----------|-----------|---|
| ICC      |           | 1/200 - 1/5000.   |
| ELISA    |           | 1/10000 - 1/100000. (Primary)   |
| WB       |           | 1/5000 - 1/50000. By chemiluminescent detection. Use at 1/5000-1/30000 with colorimetric detection. |

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors