

Anti-Human IgG4 pFc' [HP6023] ab99814

製品の概要

製品名	Anti-Human IgG4 pFc' [HP6023]
由来種	Mouse
ターゲット生物種	Human
特異性	ab99814 reacts with the pFc' portion of the heavy chain of Human IgG4.
アプリケーション	適用あり: ELISA, WB, Dot blot, IHC-Fr, ICC/IF
免疫原	Full length protein corresponding to IgG4.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	pH: 8.20 Constituent: 100% Borate buffered saline No preservatives or amine-containing buffer salts added.
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	HP6023
アイソタイプ	IgG3
軽鎖の種類	kappa

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab99814の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent dilution.

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution.
Dot blot		Use at an assay dependent dilution. (for Dot and slot immunoblotting)
IHC-Fr		Use at an assay dependent dilution.
ICC/IF		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors