

Mouse monoclonal [4E3] Anti-Human IgG1 hinge heavy chain (PE) ab99776

2 References

製品の概要

製品名	Mouse monoclonal [4E3] Anti-Human IgG1 hinge heavy chain (PE)
由来種	Mouse
ターゲット生物種	Human
特異性	ab99776 reacts with the hinge region of the heavy chain of Human IgG1 as demonstrated by ELISA. May also react with IgG from other species.
アプリケーション	適用あり: Flow Cyt, ELISA, ICC/IF
免疫原	Other Immunogen Type corresponding to IgG1.
標識	PE. Ex: 488nm, Em: 575nm

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Store In the Dark.
バッファー	pH: 7.3 Preservative: 0.1% Sodium azide Constituents: Sucrose, PBS
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	4E3
アイソタイプ	IgG1

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab99776の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		Use at an assay dependent dilution.
ELISA		Use at an assay dependent dilution.
ICC/IF		Use 0.1µg for 10 ⁶ cells.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors