

# Mouse monoclonal [hlgE2] Anti-Human IgE H&L ab435

## 製品の概要

製品名	Mouse monoclonal [hlgE2] Anti-Human IgE H&L
由来種	Mouse
ターゲット生物種	Human
アプリケーション	<b>適用あり:</b> RIA, Sandwich ELISA

## 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 6.00 Preservative: 0.095% Sodium azide Constituents: 1.47% Sodium citrate, 0.9% Sodium chloride
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	hlgE2
アイソタイプ	IgG1
軽鎖の種類	unknown

## アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab435の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
RIA		Use at an assay dependent dilution.
Sandwich ELISA		Use at an assay dependent dilution. Can be paired for Sandwich ELISA with <b>Mouse monoclonal [hlgE1] Anti-Human IgE H&amp;L (ab434)</b> . Use as Detection antibody with recommended pair.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors