

## Product datasheet

# Mouse monoclonal [JDC-12] Anti-Human lambda light chain (HRP) ab99811

★★★★★ [1 Abreviews](#) [2 References](#)

### 製品の概要

---

製品名	Mouse monoclonal [JDC-12] Anti-Human lambda light chain (HRP)
由来種	Mouse
ターゲット生物種	Human
特異性	ab99811 reacts with lambda light chains as demonstrated by ELISA.
アプリケーション	<b>適用あり:</b> ELISA, Dot blot, ICC, WB, IHC-Fr
標識	HRP

### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.40 Constituents: 50% Glycerol, 50% PBS
精製度	Affinity purified
ポリ/モノ	モノクローナル
クローン名	JDC-12
アイソタイプ	IgG1
軽鎖の種類	kappa
特記事項	

### アプリケーション

---

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab99811の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA	★★★★★ (1)	1/1000 - 1/2000.
Dot blot		Use at an assay dependent dilution.
ICC		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution.
IHC-Fr		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors