

Mouse monoclonal [KT98] Anti-Rat IgG2b H&L (HRP) ab106750

2 References [画像数 2](#)

医薬用外劇物

製品の概要

製品名	Mouse monoclonal [KT98] Anti-Rat IgG2b H&L (HRP)
由来種	Mouse
ターゲット生物種	Rat
特異性	ab106750 does not cross-react with Rat IgG1, IgG2a, IgG2c and IgM.
アプリケーション	適用あり: WB, ELISA
免疫原	Rat IgG2b
標識	HRP

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.01% Thimerosal (merthiolate) Constituent: PBS
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	KT98
アイソタイプ	IgG1

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab106750の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/1000.
ELISA		1/3000.

画像



Western blot - Mouse monoclonal [KT98] Anti-Rat IgG2b H&L (HRP) (ab106750)


All lanes : Mouse monoclonal [KT98] Anti-Rat IgG2b H&L (HRP) (ab106750) at 1/1000 dilution

Lane 1 : Rat IgG1

Lane 2 : Rat IgG2a

Lane 3 : Rat IgG2b

Lane 4 : Rat IgG2c



ELISA - Mouse monoclonal [KT98] Anti-Rat IgG2b H&L (HRP) (ab106750)

Rat immunoglobulin isotyping by ELISA using Rat primary antibodies and ab106750

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors