

Mouse monoclonal [M-4] Anti-Chicken IgM mu chain ab99715

製品の概要

製品名	Mouse monoclonal [M-4] Anti-Chicken IgM mu chain
由来種	Mouse
ターゲット生物種	Chicken
特異性	Chicken IgM (μ heavy chain specific)
アプリケーション	適用あり: ICC/IF, Flow Cyt, IP, WB

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	pH: 8.2
精製度	Purified IgM
ポリ/モノ	モノクローナル
クローン名	M-4
アイソタイプ	IgM
軽鎖の種類	kappa
特記事項	In ovo suppression of chicken B lymphocytes. In vitro activation of chicken peripheral B lymphocytes

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab99715の使用に適用されず**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		Use at an assay dependent concentration.

アプリケーション	Abreviews	特記事項
Flow Cyt		Use 1 µg for 10 ⁶ cells. ab18400 - Mouse monoclonal IgM, is suitable for use as an isotype control with this antibody.
IP		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors