

Mouse monoclonal [UHB] Anti-Human IgM mu chain (Biotin)
ab99745

製品の概要

製品名	Mouse monoclonal [UHB] Anti-Human IgM mu chain (Biotin)
由来種	Mouse
ターゲット生物種	Human
特異性	ab99745 reacts with the heavy chain of Human IgM as demonstrated by ELISA. It may also react with IgM from other species.
アプリケーション	適用あり: ICC, ELISA, ICC/IF, IHC-Fr, Dot blot, WB
標識	Biotin

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle. Store In the Dark.
バッファー	pH: 7.4 Preservative: 0.1% Sodium azide Constituent: PBS
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	UHB
アイソタイプ	IgG3

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab99745の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC		Use at an assay dependent dilution.

アプリケーション	Abreviews	特記事項
ELISA		1/1000 - 1/2000.
ICC/IF		Use 1µg for 10 ⁶ cells. for indirect immunofluorescent staining of IgM+ Human B lymphocytes.
IHC-Fr		Use at an assay dependent dilution.
Dot blot		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors