

### Rabbit Anti-Armenian hamster IgG H&L (HRP) ab5745

★★★★★ **1 Abreviews**   **10 References**   画像数 1

#### 製品の概要

製品名	Rabbit Anti-Armenian hamster IgG H&L (HRP)
由来種	Rabbit
ターゲット生物種	Armenian hamster
特異性	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase, anti-Rabbit Serum, Armenian Hamster IgG and Armenian Hamster Serum.
アプリケーション	<b>適用あり:</b> Dot blot, WB, ELISA, Electron Microscopy, IHC-P, ICC
免疫原	Armenian Hamster IgG whole molecule.
標識	HRP

#### 製品の特性


製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.50 Preservative: 0.01% Gentamicin sulphate Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA
特記事項 (精製)	BSA is IgG and protease free This product was prepared from monospecific antiserum by immunoaffinity chromatography using Armenian Hamster IgG coupled to agarose.
特記事項 (標識)	Horseradish Peroxidase (HRP).
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

**The Abpromise guarantee**   **Abpromise保証は、次のテスト済みアプリケーションにおけるab5745の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Dot blot		Use at an assay dependent concentration.
WB		1/1000 - 1/5000.
ELISA		1/25000. This product has been assayed against 1.0 ug of Armenian Hamster IgG in a standard capture ELISA using ABTS (2,2'-azino-bis-[3-ethylbenthiazoline-6-sulfonic acid]) as a substrate for 30 minutes at room temperature.
Electron Microscopy		Use at an assay dependent concentration.
IHC-P	★★★★★ (1)	1/1000 - 1/5000.
ICC		1/1000 - 1/5000.

## 画像



ab5745 tested against purified Armenian Hamster IgG. Each well was coated in duplicate with 1.0 µg of Armenian Hamster IgG. The starting dilution of antibody was 5µg/ml and the X-axis represents the log10 of a 3-fold dilution. This titration is a 4-parameter curve fit where the IC50 is defined as the titer of the antibody. Assay performed using 3% fish gelatin as blocking buffer and TMB substrate

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors