

Rabbit Anti-Donkey IgG H&L (Alkaline Phosphatase) ab6766

製品の概要

製品名	Rabbit Anti-Donkey IgG H&L (Alkaline Phosphatase)
由来種	Rabbit
ターゲット生物種	Donkey
アプリケーション	適用あり: IHC-P, WB, ELISA, Immunomicroscopy, Dot blot, ICC/IF, IHC-Fr
免疫原	donkey IgG whole molecule
標識	Alkaline Phosphatase

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.1% Sodium azide Constituents: 0.878% Sodium chloride, 1% BSA, 0.788% Tris HCl, 0.0014% Zinc chloride, 0.0095% Magnesium chloride, 50% Glycerol (glycerin, glycerine)
精製度	Affinity purified
特記事項 (精製)	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons).
特記事項 (標識)	This product was prepared from monospecific antiserum by immunoaffinity chromatography using donkey IgG coupled to agarose beads
ポリ/モノ	ポリクローナル

アプリケーション

The Abpromise guarantee **Abpromise保証は、**次のテスト済みアプリケーションにおけるab6766の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		
WB		

アプリケーション	Abreviews	特記事項
ELISA		
Immunomicroscopy		
Dot blot		
ICC/IF		
IHC-Fr		

追加情報

suitable for immunoblotting (western or dot blot), ELISA and immunohistochemistry as well as other phosphatase-antibody based enzymatic assays requiring lot-to-lot consistency

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors