

# Rabbit F(ab')<sub>2</sub> Anti-Goat IgG - H&L (Biotin), pre-adsorbed ab5753

## 1 References

### 製品の概要

|          |  |
|----------|--|
| 製品名      | Rabbit F(ab') <sub>2</sub> Anti-Goat IgG - H&L (Biotin), pre-adsorbed  |
| 由来種      | Rabbit   |
| ターゲット生物種 | Goat   |
| 特異性      | Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-biotin, anti-Rabbit Serum, Goat IgG and Goat Serum. No reaction was observed against Human Serum Proteins. |
| アプリケーション | <b>適用あり:</b> WB, ICC/IF, Dot blot, ELISA, IHC, Immunomicroscopy, Flow Cyt  |
| 吸着処理血清   | Human <a href="#">more details</a>   |
| 免疫原      | Goat IgG whole molecule  |
| 標識       | Biotin   |

### 製品の特性

|           |   |
|-----------|---|
| 製品の状態     | Liquid  |
| 保存方法      | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.  |
| バッファー     | Preservative: 0.05% Sodium azide<br>Constituents: 0.42% Potassium phosphate, 0.87% Sodium chloride, 1% BSA<br><br>BSA is IgG and protease free  |
| 特記事項 (精製) | This product was prepared from monospecific antiserum by immunoaffinity chromatography using Goat IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. |
| 特記事項 (標識) | Biotinamidocaproate N-Hydroxysuccinimide Ester (BAC) Biotin/Protein Ratio: 10-20 BAC molecules per Rabbit IgG molecule  |
| ポリモノ      | ポリクローナル   |
| アイソタイプ    | IgG   |

**The Abpromise guarantee** **Abpromise保証は、次のテスト済みアプリケーションにおけるab5753の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション         | Abreviews | 特記事項                                |
|------------------|-----------|-------------------------------------|
| WB               |           | 1/2000 - 1/10000.                   |
| ICC/IF           |           | 1/1000 - 1/5000.                    |
| Dot blot         |           | 1/10000.                            |
| ELISA            |           | 1/20000 - 1/100000.                 |
| IHC              |           | 1/500 - 1/5000.                     |
| Immunomicroscopy |           | Use at an assay dependent dilution. |
| Flow Cyt         |           | 1/200 - 1/1000.                     |

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors