

# Rabbit Anti-Human IgA H&L (HRP) ab8510

## 3 References

### 製品の概要

製品名	Rabbit Anti-Human IgA H&L (HRP)
由来種	Rabbit
ターゲット生物種	Human
特異性	This product does have a minor cross-reactivity to mouse samples.
アプリケーション	<b>適用あり:</b> Dot blot, ELISA
標識	HRP

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C.
バッファー	pH: 7.50 Constituents: 0.164% Sodium phosphate, 50% Glycerol (glycerin, glycerine), 0.58% Sodium chloride
精製度	Immunogen affinity purified
特記事項 (精製)	From the crude polyclonal the crossreactive antibodies were extracted by incubation with Sepharose-bound human IgM and IgG. Specific antibodies were absorbed by incubation with Sepharose-bound human IgA. Specific antibodies were eluted by acidic buffer at pH 2.5 followed by neutralisation and dialysis. After repeated binding with immobilized human IgA a minimum of 65% protein bound.
特記事項 (標識)	The purified polyclonal was conjugated to horseradish peroxidase according to the periodate method followed by gel-filtration and ion-exchange chromatography to clear unbound conjugate.
ポリモノ	ポリクローナル

### アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab8510の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Dot blot		1/100.
ELISA		1/1000.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors