

Rabbit monoclonal [MlgG-R180-9-8] Anti-Mouse IgG1+IgG2 H&L ab125916

リコンビナント **RabMAb**

画像数 1

製品の概要

| | |
|----------|---|
| 製品名 | Rabbit monoclonal [MlgG-R180-9-8] Anti-Mouse IgG1+IgG2 H&L |
| 由来種 | Rabbit |
| ターゲット生物種 | Mouse |
| 特異性 | ab125916 is specific for Mouse IgG1 and Mouse IgG2a. The antibody does not recognise Mouse IgG2b, Mouse IgG3, Mouse IgA, Mouse IgM, Human IgG or Rat IgG. |
| アプリケーション | 適用あり: ELISA |
| 免疫原 | Other Immunogen Type corresponding to Mouse IgG1 + IgG2a. |

製品の特性

| | |
|--------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| バッファー | pH: 7.20 Preservative: 0.09% Sodium azide Constituent: 99% PBS |
| 精製度 | Protein A purified |
| ポリ/モノ | モノクローナル |
| クローン名 | MlgG-R180-9-8 |
| アイソタイプ | IgG |
| 特記事項 | <p>ab125916 has weak cross-reactivity with mouse IgG1.</p> <p>Human, Rat: We have preliminary internal testing data to indicate this antibody may not react with these species. Please contact us for more information.</p> <p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply |

- Animal-free production

For more information [see here](#).


Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to [RabMAb[®] patents](#).

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab125916の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|---------|
| ELISA | | 1/1000. |

画像



Direct ELISA antibody dose-response curve using ab125916 between 0-1000 ng/ml. Antigen concentration of 1000 ng/mL. An alkaline phosphatase-conjugated goat anti-rabbit IgG (H+L) (1/2500) was used as the secondary antibody.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors