

Rabbit Anti-Mouse IgG H&L ab6709

50 References 画像数 1

製品の概要

製品名	Rabbit Anti-Mouse IgG H&L
由来種	Rabbit
ターゲット生物種	Mouse
アプリケーション	適用あり: WB, IP, Immunodiffusion, Conjugation, ICC/IF, ELISA, IHC-P, IHC-Fr
免疫原	Mouse IgG, whole molecule

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.01% Sodium azide Constituents: 0.42% Potassium phosphate, 0.87% Sodium chloride
精製度	Affinity purified
特記事項 (精製)	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Mouse IgG coupled to agarose beads
ポリ/モノ	ポリクローナル


アプリケーション

The Abpromise guarantee      Abpromise保証は、次のテスト済みアプリケーションにおけるab6709の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/2000 - 1/10000.
IP		Use at an assay dependent concentration.
Immunodiffusion		Use at an assay dependent concentration.

アプリケーション	Abreviews	特記事項
Conjugation		Use at an assay dependent concentration.
ICC/IF		1/1000 - 1/5000.
ELISA		1/20000 - 1/400000.
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.

## 画像



Sensitivity of ab6709 against mouse IgG.

The starting dilution of ab6709 was 5µg/ml and the X-axis represents the Log10 of a 3-fold dilution. This titration is a 4-parameter curve fit where the IC50 is defined as the titer of the antibody. The blocking agent was 3% fish gelatin. The detection antibody was a Goat anti-Rabbit IgG Antibody Peroxidase conjugate. TMB color development.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors