

# Rabbit Anti-Mouse IgM mu chain ab9175

## 1 References

### 製品の概要

製品名	Rabbit Anti-Mouse IgM mu chain
由来種	Rabbit
ターゲット生物種	Mouse
特異性	The affinity purified anti IgM mu antibodies only react with IgM by immunodiffusion and Immunoelectrophoretic techniques. A 1 mg/ml solution of the antibody is reacted against normal mouse serum and IgG1, IgG2a, IgG2b, IgG3, IgM and IgA. By ELISA, there is minimal reactivity to IgG.
アプリケーション	<b>適用あり:</b> ELISA, Immunoelectrophoresis, WB
免疫原	Rabbits were immunized with highly purified normal mouse IgM in Freund's adjuvant.

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.20
精製度	Immunogen affinity purified
特記事項 (精製)	The antiserum was solid phase absorbed to make it IgM mu chain specific and the anti mu chain antibodies were isolated by solid phase immunoaffinity chromatography.
ポリ/モノ	ポリクローナル

### アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab9175の使用に適用されます**  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.
Immunoelectrophoresis		Use at an assay dependent concentration.

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent concentration.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors