

# Rabbit Anti-Sheep IgG H&L (TRITC) ab6744

## 1 References

### 製品の概要

製品名	Rabbit Anti-Sheep IgG H&L (TRITC)
由来種	Rabbit
ターゲット生物種	Sheep
アプリケーション	<b>適用あり:</b> Immunomicroscopy, Flow Cyt, ICC/IF, IHC-P, IHC-Fr, ELISA
免疫原	Sheep IgG whole molecule
標識	TRITC. Ex: 547nm, Em: 572nm

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.878% Sodium chloride, 1% BSA, 0.424% Tripotassium orthophosphate
精製度	Affinity purified
特記事項 (精製)	.This product was prepared from monospecific antiserum by immunoaffinity chromatography using Sheep IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
特記事項 (標識)	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Sheep IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities
ポリモノ	ポリクローナル
特記事項	Tetramethylrhodamine isothiocyanate (TRITC) (Molecular Weight 444 daltons) Absorption Wavelength: 550 nm Emission Wavelength: 570 nm Fluorochrome/Protein Ratio: 2.8 moles TRITC per mole of Rabbit IgG.

### アプリケーション

The Abpromise guarantee

**Abpromise保証は、** 次のテスト済みアプリケーションにおけるab6744の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
Immunomicroscopy		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent dilution.
ICC/IF		1/1000 - 1/5000.
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
ELISA		1/10000 - 1/50000.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors