



Product datasheet

Rabbit Anti-Sheep IgM H&L (HRP) preadsorbed ab112763

2 References

製品の概要

| | |
|----------|--|
| 製品名 | Rabbit Anti-Sheep IgM H&L (HRP) preadsorbed |
| 由来種 | Rabbit |
| ターゲット生物種 | Sheep |
| 特異性 | By immunoelectrophoresis and ELISA ab112763 reacts specifically with sheep IgM. May cross react with IgM from other species. |
| アプリケーション | 適用あり: WB, ELISA, IHC-P, ICC/IF |
| 免疫原 | Other Immunogen Type corresponding to IgM. |
| 標識 | HRP |

製品の特性

| | |
|-----------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C. |
| バッファー | pH: 6.8 Constituents: 0.2% BSA, 0.05% CMIT/MIT based preservative |
| 精製度 | Immunogen affinity purified |
| 特記事項 (精製) | Antiserum was solid phase adsorbed to ensure class specificity. ab112763 was isolated by affinity chromatography using antigen coupled to agarose beads. |
| 特記事項 (標識) | Molar enzyme/antibody protein ratio is 4:1. |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab112763の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|------|
| | | |

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|--|
| WB | | 1/1000 - 1/10000. colorimetric detection; 1/1000 - 1/30000 chemiluminescent detection. |
| ELISA | | 1/10000 - 1/100000. as a primary. |
| IHC-P | | 1/200 - 1/500. |
| ICC/IF | | 1/200 - 1/500. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors