

# Rabbit Anti-Syrian Hamster IgG H&L (HRP) ab6783

## 1 References

### 製品の概要

製品名	Rabbit Anti-Syrian Hamster IgG H&L (HRP)
由来種	Rabbit
ターゲット生物種	Syrian hamster
アプリケーション	<b>適用あり:</b> ICC, Dot blot, IHC-P, IHC-Fr, Immunomicroscopy, WB, ICC/IF, ELISA
免疫原	hamster IgG whole molecule
標識	HRP

### 製品の特性

製品の状態	Lyophilized:Store vial at 4° C prior to restoration. Restore with 1.0 mL of deionized water (or equivalent). For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. Expiration date is one (1) year from date of restoration.
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.01% Gentamicin sulphate Constituents: 0.87% Sodium chloride, 1% BSA, 0.42% Potassium phosphate
精製度	IgG fraction
特記事項 (精製)	Horseradish Peroxidase (HRP).
特記事項 (標識)	This product was prepared from monospecific antiserum by immunoaffinity chromatography using hamster IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities
ポリ/モノ	ポリクローナル

### アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、**次のテスト済みアプリケーションにおけるab6783の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC		Use at an assay dependent dilution.
Dot blot		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
Immunomicroscopy		Use at an assay dependent concentration.
WB		1/1000 - 1/5000.
ICC/IF		1/500 - 1/2500.
ELISA		1/16000 - 1/64000.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors